
	DEPARTMENT PARASITOLOGY BPRC Rijswijk	SOP: 001 Version: 001
	Enzyme Linked ImmunoSorbent Assay (ELISA)	Page 1 of 3 Date: 17/11/10

Index	Page
A.1. Aim	2
B.1. Equipment and instruments	2
B.2. Materials	2
C.1. Procedure	2
D.1. Schematics	3
D.2. Example on pre-dilutions	3

	DEPARTMENT PARASITOLOGY BPRC Rijswijk	SOP: 001 Version: 001
	Enzyme Linked ImmunoSorbent Assay (ELISA)	Page 2 of 3
		Date: 17/11/10

A.1. Aim

The detection and quantification of antibodies (Ab) by a sandwich Enzyme Linked ImmunoSorbent Assay. This assay can measure Ab titers from 0.78 ng/ml to 100 ng/ml from serum, supernatant or plasma, depending on which coating is used.

B.1. Equipment and instruments


- ELISA plates (Greiner microlon #655092)
- 50 ml tubes: Greiner Bio-one #22761
- Pipettes: 2 ml: ALP #PN2E1
5 ml: ALP #PN5E25
10 ml: ALP #PN10E25
25 ml: ALP #PN25E1
- Deep well 96-well plates 1ml (NUNC, #278605)
- Plate washer (Bio-Tek instruments, ELx405)
- Plate reader (BIORAD microplate reader, model 680)

B.2. Materials

- Coating antigen : AMA-1, MSP-1, DICO 1,2,3 or Goat anti Hu, Rb, Ra or Mo.
- Coating Buffer : PBS (Gibco #10010-031)
- Blocking Buffer : PBS + 0,05% v/v Tween-20 + 3% w/v BSA
Tween-20: Sigma # P1379
BSA: Sigma #A9647
- Wash Buffer : PBS + 0,05% v/v Tween-20
- Dilution Buffer : PBS + 0,05% v/v Tween-20 + 0,5% w/v BSA
- Standard e.g. pool of final bleed of animals immunized with AMA-1, MSP-1, DICO 1,2,3 or total Hu, Rb, Rt or Mo IgG.
- Anti Hu, Rb, Ra or Mo IgG or IgM conjugated to Alkaline Phosphatase
- DEA buffer
Diethanolamine: Sigma #D83303
MgCl₂*6H₂O: Sigma #M0250
To make 500 ml: 0,15% w/v MgCl₂*6H₂O in 500 ml MQ (e.g. 0,75 g) + 492 µl Diethanolamine pH 9,8
- Substrate PNPP (Para Nitro Phenyl Phosphate hexahydrate) (Fluka cat# 71768)

C.1. Procedure

1. Coat ELISA plates with the desired coating, see schedule below for concentrations.
Make dilutions in Coating Buffer and add 100 µl per well.
2. Incubate overnight (O/N) at 4°C covered with a plastic seal or lid.
3. Remove coating from plates by inverting the plates with a vigorous wrist action.
4. Block with 200 µl/well of Blocking Buffer.
5. Incubate 1 hour at Room Temperature.
6. Wash with plate washer program 9 (each well aspirated, filled with wash buffer and the wash buffer is aspirated out again. This cycle is performed five times).
7. Prepare samples and standards:
Dilute pre immune samples 10 - 100 times in Dilution buffer as starting dilution and final bleed samples 1000 – 10000 times.
Dilute the standard and make a curve starting at an IgG concentration of approximately 50 ng/ml.
Example on how to make predilutions, check the end of the protocol
8. Add 100 µl of Dilution Buffer to all wells.

	DEPARTMENT PARASITOLOGY BPRC Rijswijk	SOP: 001 Version: 001
	Enzyme Linked ImmunoSorbent Assay (ELISA)	Page 3 of 3
		Date: 17/11/10

9. Add 100 µl of each pre diluted sample or standard IgG sample to the first wells of a column and make a 3-fold serial dilution over 8 wells using 50 µl. The standard IgG sample is serial diluted in 2-fold over 7 wells, using the last well as a blanc.
10. Incubate 1 hour at RT.
11. Wash with plate washer program 9.
12. Add 100 µl conjugate (Alkaline Phosphatase) in the right dilution (see schedule below).
13. Incubate 1 hour at RT.
14. Wash with plate washer program 9.
15. Add 100 µl of 1 mg/ml PNPP in DEA buffer.
16. Incubate 30 minutes at RT.
17. Read OD at 405 nm on a Microplate Reader.

D.1. Schematic:

Action	Material 1	Supplier	Concentration	Dilute in	Incubate
Coating	AMA-1, MSP-1	BPRC	1 µg/ml	Coating buffer = PBS	O/N 4 C
	Dico 1, 2 or 3				
	Sh anti Rb IgG	Sigma R-3631	1 : 4000		
	Go anti Mo IgG	Thermo S. 31160	1 : 8000		
	Go anti Ra IgG	Pierce 31220	1 : 2000		
	Go anti Rb IgG	Thermo S. 31210	1 : 8000		
	Go anti Hu IgG	Thermo S. 31130	1 : 8000		
	Go Anti Hu IgM	Sigma I-2386	5 ug/mL		
Blocking	PBS-Tween-BSA	Merck	0,05% Tween	PBS	1 hr RT
		8.22181.0500 Sigma A-9647	3% BSA		
IgG Standard	Human IgG	Sigma I-4506	100ng/mL – 2 fold	Dilution buffer, after reconstitution in 150mM NaCl	
	Mouse IgG	Sigma I-5381	100ng/mL – 2 fold		
	Rabbit IgG	Sigma I-5006	100ng/mL – 2 fold		
Samples	Pre immune sera		Start 10 – 100x	Dilution buffer	
	Immune sera		Start 1000 – 10000x		
Conjugate	Go anti Rb IgG-AP	Pierce 31340	1:1250	Dilution buffer	1 hr RT
	Go anti Ra IgG-AP	Pierce 31350	1:1250		
	Go anti Mo IgG-AP		1:1250		
	Go anti Hu IgG-AP	Pierce 31310	1:1250		
	Go anti Hu IgM-AP	Sigma A-2189	1: 10000		
Substrate	PnPP	Fluka 71768	1 mg/ml	DEA buffer	30 min RT

D.2. Example on how to make pre-dilutions:

If the first dilution in the well should be 50ng/mL followed by 2 fold dilution series over seven wells, then the pre-dilution needs to be twice more concentrated, so 100ng/mL. After adding 50ul dilution buffer to each well, add 50ul of the pre-dilution (100ng/ml) to the top wells. Which dilutes it to 50ng/ml. Then 50ul is taken out to be resuspended in the next well, etc, to make the dilution series.