



Enzyme-Linked Immunosorbent Assay (ELISA)

1.0 REAGENTS

1.1. PHOSPHATE-BUFFERED SALINE (PBS)

2mMKH₂PO₄

4mM Na₂HPO₄

0.15M NaCl

pH to 7.5

1.2. ANTIGEN DILUTION BUFFER

0.1 M carbonate buffer (Na₂CO₃/NaHCO₃)

1.3. COATING BUFFER

15mM Na₂CO₃

35mM NaHCO₃

pH 9.3

1.4. WASHING BUFFER PBS-TWEEN-20

PBS as above +0.05% Tween-20

1.5. BLOCKING BUFFER

1% [w/v] skimmed milk powder in washing buffer

1.6. HORSERADISH PEROXIDASE-CONJUGATED RABBIT SERUM

Dako Ltd, High Wycombe, UK

1.7. SUBSTRATE

(0.04mg ml⁻¹ o-phenylenediamine [Sigma prod no. 78410] + 0.012% H₂O₂)

1.8. COATING BUFFER

15mM Na₂CO₃

35mM NaHCO₃
pH 9.3

1.9. STOP SOLUTION

2M H₂SO₄

1.10. DEVELOPING BUFFER

24.5mM citric acid monohydrate and 52mM Na₂HPO₄, pH 5.0).

2.0 EQUIPMENT

2.1. 96-WELL MICRITITRE PLATES

Immulon 4 HBX, Dynex

2.2. ELISA PLATE WASHER

Ultrawash Plus, Dynex

2.3. ELISA READER

3.0 METHOD

1. Design the layout of the experiment in the 96-well format. All test sera should appear twice on the same plate. Included on each plate should be duplicate wells containing a) reagents alone and b) known negative sera controls
2. Add 100 µL of antigen made up to the desired concentration e.g. 0.5 µg/mL in 0.1 M carbonate buffer (Na₂CO₃/NaHCO₃) per to 96-well microtitre plates (Immulon-4; Dynatech, Billingshurst, UK) and incubate at room temperature overnight.
3. Remove antigen and wash plates 3x with Washing buffer using the Ultrawash Plus Plate washer. Shake dry and add 200 µL of blocking buffer (1% milk powder in PBS, 0.05% Tween 20) and incubate at room temperature for 5 hours.
4. Wash plates 3x with Washing buffer and add 100 µL of serum (diluted 1/1000 in PBS) each of duplicate wells and incubate overnight at 4°C.
5. Wash plates, shake dry and add 100 µL of a 1/6000 dilution of horseradish peroxidase conjugated rabbit anti-human IgG antibody (Dako Ltd) and incubate for 3 h at room temperature. Develop the reaction with 200 µL of H₂O₂ and o-phenylenediamine at 4°C for 10 min. Stop the reaction by adding 25 µL of 2 M H₂SO₄ per well, and read the

optical density (OD) was measured at 492 nm using aELISA reader. For each tested serum, the OD